



EGG INCUBATOR

Outer wall is made up of M.S. sheet duly painted and inside of S.S. . The gap between inner and outer wall is filled with high grade mineral glass wool. Temp. range from ambient to 70° C is controlled by “JUMO”/ EGO” German capillary type thermostat with an accuracy of + 0.5° c. The unit is provided with manual tilting device, so as to tilt the tray holding the eggs. To work on 220/230 volts AC supply

CAPACITY

- (a) 100 eggs
- (b) 150eggs
- (c) 200eggs
- (d) 250eggs
- (e) 600eggs

Manufacturers & Suppliers:

LABPRO INTERNATIONAL

48 A 1ST FLOOR, KUNJ VIHAR, BEHIND B.D. FLOOR MILL, AMBALA CANTT-133001 HARYANA (INDIA)



PH. 0171-2690847 (off.), MOB. 9034101997(K), MOB. 9728790667(K)

E-mail: labprointernational@gmail.com , labpro@india.com