



Heating Mantle

We offer an assortment of heating mantles that have gained huge appreciation in the market for their host of qualitative features. These instruments are manufactured using optimum quality raw material sourced from reliable vendors.

Technical Specifications:

- Body Made from Aluminium Sheet (Epoxy Coated)
- Mantle Made From Glass Yarn
- Sizes: 250ml, 500ml, 1000ml, 2000ml, 3000ml, 5000ml, 10000ml
- Temp. up-to 350 °c Energy Regulator
- Works on 220 V AC