



### SSU-133 V.D.R.L. ROTATOR

Used for V.D.R.L. Tests, blood grouping tests, various agglutination tests or for mixing of solutions in small bottles, flasks or beakers. It has a platform of size 300 x 300 mm. and can accommodate blood bottles, slides or flasks using a spring holder. Supplied complete with cord and plug to work on 220/230 volts A.C.,

- (a) Fixed speed (fixed at 180 RPM)
  - (b) Variable speed
- OPTIONAL ACCESSORIES
- (a) Automatic 0-60 minutes "GERMAN" Timer

## LABPRO INTERNATIONAL

48 A 1<sup>ST</sup> FLOOR, KUNJ VIHAR, BEHIND B.D. FLOOR MILL, AMBALA  
CANTT-133001 HARYANA (INDIA) PH. 0171-2690847 (off.), MOB.  
9034101997(K), MOB. 9728790667(K)

E-mail: [labprointernational@gmail.com](mailto:labprointernational@gmail.com) , [labpro@india.com](mailto:labpro@india.com)

*Manufacturers & Suppliers:* LABORATORY EQUIPMENTS, GLASSWARE ITEMS, PHARMACY, CHEMISTRY EQUIPMENTS, SCIENTIFIC EQUIPMENTS, SCHOOL / COLLEGE / UNIVERSITY EQUIPMENTS, INDUSTRIES EQUIPMENTS etc