



DISTILLATION UNIT

Available in output capacity of 4ltr/hr only it is very rugged and has a very long life as compared to the routine water still. All parts such as boiler, top cover and JOINTLESS condenser can be frequently cleaned by opening the screws of removable condenser. To work on 220/230 volts A.C.

CAPACITY LOAD

With cutoff system

LABPRO INTERNATIONAL

**48 A 1ST FLOOR, KUNJ VIHAR, BEHIND B.D. FLOOR MILL, AMBALA
CANTT-133001 HARYANA (INDIA)**

PH. 0171-2690847 (off.), MOB. 9034101997(K), MOB. 9728790667(K)

E-mail: labprointernational@gmail.com , labpro@india.com

Manufacturers & Suppliers:

LABORATORY EQUIPMENTS,